


Abstract 6929: Inter- and intratumoral *PIK3CA* subclonal diversity in breast cancer contextualized by single-cell multiomics

J. Zawistowski¹, I. Salas-Gonzalez¹, T. Tate¹, T. Morozova¹, K. Kennedy¹, D. Arvapalli¹, J. Remington¹, J. Marks², E. Hwang², G. Harton¹, V. Weigman¹, J. A. A. West¹

¹BioSkrby Genomics, Durham, NC, United States, ²Duke University Medical Center, Durham, NC, United States


Introduction




The diversity of molecular determinants of the transition from the premalignant state of DCIS to invasive cancer are multi-faceted and will require single-cell, multimodal analysis across the Central Dogma of Molecular Biology

- 55,720 new cases of DCIS in 2023¹
 - 20-50% of DCIS progress to invasive breast cancer²
 - Virtually all patients with DCIS have surgery; almost a third of patients have unilateral or bilateral mastectomy³
 - Identifying markers of progression could reduce overtreatment of low-risk lesions
- Is it possible to "Treat the biology; not the diagnosis"?


Methods



ResolveOME™ single-cell workflow to obtain unified genomic, transcriptomic, and targeted protein data.
A. Sample singulation workflow upstream of ResolveOME™ reactions. **B.** ResolveOME™ workflow steps and time to completion (top) and performance metrics for each analyte type (bottom). **C.** Inclusion of BioLegend oligo-conjugated antibodies in ResolveOME™ workflow to report surface protein expression in conjunction with genomic/transcriptomic tiers.



DNA Performance Characteristics	Observed Values	RNA Performance Characteristics	Observed Values	Protein Performance Characteristics	Observed Values
Accuracy	99.99%	Protein coding genes detected	3,451 ± 731.71*	Surface proteins detected	~ 165
Sensitivity	96.04%	Proportion exonic	0.77	Antibody-derived tag (barcode) signal (1M reads)	10 ³ - 10 ⁴
Specificity	99.99%	Average concordance	0.97	Isotype control signal (1M reads)	< 10 ²
Precision	98.75%	Reproducibility (CV)	32.9%		
Allelic balance	93.46%				
Genome coverage	97.59%				



Results




Figure 1: *PIK3CA* subclonal diversity in DCIS/IDC tumors. **A.** cBioPortal lollipop diagram showing domain location of *PIK3CA* mutations (driver, or variants of unknown significance, VUS) among invasive breast cancer samples of published studies; lollipop height indicates mutation frequency among cohorts. Highlighted are domain locations of mutations found in single cells from Patient A ($\Delta E109$, K111E, H1047R) and Patient B (N345T) of our study.







Figure 2: Comprehensive single-cell CNV profiling of 12 DCIS/IDC tumors. Estimated ploidy for each chromosome is indicated for single cells (rows) from all patient tumors of the study, including Patient A and B harboring *PIK3CA* missense mutations. Cell typing ID is shown for each cell as well as pathological grouping of each cell.

Patient #	Cell #	<i>PIK3CA</i> mutation	TCGA tissue	BaseJumper cell type	HPCA cell type	Cell cycle stage
A	1	K111E	breast	epithelial	epithelial	S
A	2	H1047R	breast	epithelial	epithelial	G1
A	3	H1047R	breast	epithelial	epithelial	S
A	4	H1047R	breast	epithelial	epithelial	G1
A	5	$\Delta E109$	breast	epithelial	epithelial	G1
A	6	$\Delta E109$	breast	epithelial	epithelial	S
A	7	$\Delta E109$	thyroid	epithelial	epithelial	G2M
A	8	$\Delta E109$	breast	endothelial	endothelial	S
A	9	$\Delta E109$	breast	hepatocyte	epithelial	S
A	10	$\Delta E109$	breast	epithelial	epithelial	G1
A	11	$\Delta E109$	breast	epithelial	epithelial	S
A	12	$\Delta E109$	breast	epithelial	epithelial	G1
A	13	$\Delta E109$	breast	epithelial	epithelial	G2M
A	14	$\Delta E109$	breast	epithelial	epithelial	G1
A	15	$\Delta E109$	breast	epithelial	epithelial	G1
A	16	$\Delta E109$	breast	epithelial	epithelial	G2M
B	1	N345T	lung	macrophage	epithelial	G1
B	2	N345T	-	-	-	-



Conclusions

- A *PIK3CA* mutational spectrum, displaying inter- and intra-tumoral variation, is defined within single cells from a subset of patients in a DCIS/IDC cohort
- Distinct copy number alteration profiles are associated with specific *PIK3CA* mutations, suggestive of a cooperativity between the two in clonal evolution and in invasiveness capacity
- RNA and surface protein-level data contextualize the genomic lesions of *PIK3CA* mutations and copy number alterations by defining cell identity and cell state

References

1. Siegel et al. CA Cancer J Clin (2023) <https://doi.org/10.3322/caac.21763>
2. Cowell et al. Mol Onc (2013) <https://doi.org/10.1016%2Fj.molonc.2013.07.005>
3. Ward et al. CA Cancer J Clin (2015) <https://doi.org/10.3322/caac.21321>
4. Marks et al. bioRxiv (2023) <https://doi.org/10.1101/2022.04.29.489440>
5. Gorringer et al. Modern Pathology (2015) <https://doi.org/10.1038/modpathol.2015.75>

Acknowledgements

We thank the patients involved in this study for their gracious gift of tissue to fuel discovery: IRB PRO00034242 "Biologic Characterization of the Breast Cancer Tumor Microenvironment" (PI: Hwang).