

Enrichment of viable cells from patient samples with LeviCell technology as input into the ResolveOME™ single cell multiomic workflow

Alison Rojas, Swetha D. Velivela, Jessica N. Borja, Jesse A. Ramirez, Tatiana V. Morozova, Jeff G. Blackinton, Jon S. Zawistowski


BioSkryb Genomics, Durham, NC, USA




Introduction

Single cell multiomic analysis is at the forefront of studies driving insight into tumor clonal evolution and somatic mosaicism in normal tissues, as well as into edited cell line surveillance. The ability to generate unified genomic and transcriptomic information with fidelity is wholly dependent on the input of viable cells into the amplification chemistry, otherwise concurrent genotype and phenotype are not accurately captured. We present the coupling of the Levitas LeviCell 1.0 instrument, where viable vs dying/dead cells differentially levitate in a microfluidic cartridge within a magnetic field, to single cell dispensing technologies providing input into ResolveOME multiomic amplification technology. Primary cryopreserved bronchoalveolar lavage cells were either directly inputted into FACS or HP D100 instruments for dispensing into 96 well PCR plates or inputted into the LeviCell 1.0 prior to dispensing. For FACS dispensing, even with the use of Calcein-AM and propidium iodide dual viability staining, levitation of primary cells (2 independent samples, 2 replicates each) prior to dispensing increased the percentage of cells yielding genomic amplification product from 51.25% to 88.75% and, for transcriptomic yield, from 52.50% to 97.50%. A similar trend of reduced reaction yield dropouts was obtained for HP D100 dispensing for both genomic amplification and cDNA yield, albeit with a more modest mean improvement of 17%—yet LeviCell viability enrichment importantly reduced the amount of sporadic clogging of the HP D100 microfluidic cartridge. Low-coverage genomic sequencing metrics assessing library complexity and amplification uniformity were overall comparable between pre- and post-levitation, as were the transcriptomic sequencing metrics of mitochondrial transcript percentage, exonic:intergenic read ratio, and number of expressed genes detected. The marked reduction of reaction failures or sub-optimal reaction yields while maintaining robust sequencing metrics with this LeviCell-dispensing-ResolveOME workflow is a facile solution for maximizing the utility of rare and valuable samples, as well as for maximizing the success rate of a given ResolveOME run with often challenging clinical samples.


Introduction



Strategy employed for viability assessment of single cells with FACS (left). A dual method of positive selection with Calcein AM and negative selection with propidium iodide maximizes viable cell output from FACS. **LevitasBio LeviCell 1.0 system (right).** Suspended cells flow through microfluidic cartridge subjected to a magnetic field and a contrast agent, enriching viable cells through buoyancy properties. The bottom image shows separation of viable and non-viable murine T-cell lymphoma (EL4) cells corroborated with acridine orange and propidium iodide viability dyes.



Methods



Conclusions

- For patient samples, which often are of poor integrity, levitation prior to FACS or cell dispensing significantly improved the fraction of successful single-cell genomic and transcriptomic amplification reactions from ResolveOME multiomic chemistry
- DNA and RNA performance metrics were overall concordant between pre- and post-levitation samples
- For precious clinical samples, first enriching for viability by levitation, even in the presence of viability gating in FACS, maximizes ResolveOME reaction output

Results









Figure 3: ResolveOME genomic amplification of single cells in the presence or absence of LeviCell 1.0 viability enrichment. FACS-sorted single airway cells were subjected to the ResolveOME multiomic workflow following LeviCell 1.0 treatment (Post-LeviCell) or directly dispensed in the absence of levitation (Pre-LeviCell). Amplification yield of the genomic arm of the ResolveOME multiomic chemistry is shown for both conditions, across two operators. Levitation significantly reduced the number of zero- or low-yielding ResolveOME reactions for both operators.



References

- Marks JR et al. Unifying comprehensive genomics and transcriptomics in individual cells to illuminate oncogenic and drug resistance mechanisms. bioRxiv preprint: doi: <https://doi.org/10.1101/2022.04.29.4894402>.
- Gonzalez-Pena, V. et al. Accurate genomic variant detection in single cells with primary template-directed amplification. *Proc Natl Acad Sci U S A* 118, e2024176118 (2021).

Acknowledgements

We thank Drs. Peter Park, Lovelace Luquette, and Craig Bohrsen (Harvard Medical School) for collaboration with airway lavage samples and thank the patients involved in this study for their gracious gift of tissue to fuel discovery. We thank Katy Richards-Hrdlicka, Brian Kierce, and Joseph Breier at LevitasBio for their collaboration in generating these data.