

[Return to website >](#)

A revolution in **resolution**


NOW AVAILABLE COMMERCIALY

ResolveDNATM

Whole genome sequencing solution

BioSkryb
GENOMICS

Cells explored. Answers revealed.




Powering a genomics revolution with superior resolution

Massively parallel bulk sequencing has vastly expanded our knowledge of the living world. However, higher resolution is needed to understand and leverage the contributions of individual cells to the biology of organisms and ecosystems.

BioSkryb is revolutionizing the field of single-cell genomics with PTA (Primary Template-directed Amplification), a novel technology that enables accurate and reproducible whole genome amplification (WGA). PTA recovers >95% of the genomes of single cells with higher fidelity and uniformity than other WGA methods. This industry-leading performance is highly scalable, supporting robust, reliable and routine single-cell and ultra-low input genomics applications in clinical, translational and life sciences research.

BioSkryb is developing complete sample-to-analysis workflow solutions for variant detection and characterization from single-cell and ultra-low input samples. Our PTA-based ResolveDNA™ Whole Genome Amplification Kits form the core of whole genome and targeted sequencing workflows, enabling single-nucleotide variant (SNV) and copy number variant (CNV) analysis in cancer genomics, cardiology, neurology, immunology, toxicology and preimplantation genetic testing.

BioSkryb
G E N O M I C S



Application Areas

BioSkryb core technology supports a broad range of single-cell and ultra-low input applications



Cancer genomics



Prenatal genetic testing (PGT)



Cardiology



Microbiome research



Neurology




Toxicology



Immunology




Bioprocessing




BioSkryb
GENOMICS

ResolveDNA™ Whole Genome Sequencing Workflow

The ResolveDNA™ Whole Genome Sequencing Workflow is compatible with single cells, multiple cells and low-input (>4 pg to 10 ng) DNA samples. Whole genome amplification (WGA) is performed in three easy steps. WGA products are converted to libraries for Illumina® sequencing. The BaseJumper Bioinformatics Platform offers automated data processing and convenient visualization of variants.



Cell Sorting



The ResolveDNA™ Whole Genome Sequencing Workflow is compatible with fluorescence-activated cell sorting (FACS), as well as microfluidic and droplet-based cell sorting methods.

The ResolveDNA™ FACS Kit contains 96-well plates and sealing film, as well as ResolveDNA PTA-Grade Cell Buffer—specifically optimized for PTA-based whole genome amplification from sorted cells.



FACS Kit components are also available as part of ResolveDNA Starter Packs. The Complete Starter Pack includes additional validated reagents, consumables and equipment for the ResolveDNA Whole Genome Sequencing Workflow. A Consumables Only Starter Pack is also available.

[Product Details >](#)

BioSkryb
GENOMICS

Whole Genome Amplification

FEATURES & BENEFITS



ResolveDNA™ Whole Genome Amplification Kits overcome the inherent challenges of low coverage uniformity in single-cell genomics that lead to increased sequencing costs and complex data analysis.

The controlled reaction parameters employed in this PTA-based kit enables reproducible recovery of >95% of the genomes of single cells and limited DNA input samples with industry-leading uniformity and accuracy.

Key features and benefits of ResolveDNA Whole Genome Amplification Kits include:

- Significantly lower allelic dropout and biases compared to existing WGA methods that yield low and variable coverage across the genome
- Specific amplification of the primary template with >97% of reads mapping to the human genome and no detectable product in no template control reactions
- Simple, user-friendly workflow that requires less than 45 mins of hands-on time
- Consistent WGA fragment sizes and yields from picogram to nanogram inputs

[Product Details >](#)

BioSkryb
GENOMICS

Whole Genome Amplification

WORKFLOW



Compatible with:

- Single or multiple cells; ultra-low amounts of DNA (>4 pg—10 ng)
- Any cell collection method, plate or tube format

Three easy steps: cell lysis, whole genome amplification, bead-based cleanup of amplified DNA


WGA products are converted to libraries for **multiplexed Illumina® sequencing** (no fragmentation required)



Typical single-cell whole genome sequencing workflow. Steps highlighted in light blue are performed with the ResolveDNA Whole Genome Amplification Kit. The bead-based cleanup is performed with the Resolve DNA Bead Purification Kit.

Whole Genome Amplification

LIBRARY QC METRICS




ResolveDNA™ Whole Genome Amplification Kits offer consistent, low-microgram yields from single cell inputs.

This enables robust library prep workflows, and allows for repeat analysis from the same sample. Reproducible, low-nanogram yields from no template controls (NTCs) confirm that only the product of interest is amplified.

Whole Genome Amplification

LIBRARY QC METRICS



ResolveDNA™ Whole Genome Amplification Kits enable control over fragment size.


Unlike other WGA methods, the innovative PTA chemistry and quasi-linear isothermal process prevents amplification of daughter amplicons. This results in relatively short (~250 bp to >1.5 kb) reaction products. No fragmentation is required prior to library preparation, but bead-based size selection after library amplification may be used to obtain the optimal library insert size for your sequencing read length/application.

Whole Genome Amplification


SEQUENCING PERFORMANCE



ResolveDNA™ WGA Kit



Single cell MDA



ResolveDNA Whole Genome Amplification Kits offer superior coverage and uniformity in single cell WGS.

WGA was performed with the ResolveDNA WGA Kit (left) or single cell MDA (right). Plots show a portion of chromosome 1 (100 kb bins). The central area (poorly covered with both methods) corresponds to the centromere.

Whole Genome Amplification

SEQUENCING PERFORMANCE



Method	ResolveDNA	Mixed Method A	MDA A	MDA B	Mixed Method B	Mixed Method C	DOP-PCR
Genome mapping	97%	91%	88%	55%	88%	55%	52%
Genome recovery	97%	73%	65%	59%	50%	33%	20%
CV of coverage	0.8	1.3	1.8	2.3	2.6	3.2	3.6
SNV sensitivity	76%*	38%	40%	35%	25%	10%	4%
SNV specificity	94%	92%	93%	92%	57%	65%	44%

Sensitivity and specificity are based on positions that are 15X coverage *SNV sensitivity for ResolveDNA approaches 85% at 30X coverage.
CV: Coefficient of Variation SNV: Single Nucleotide Variation

ResolveDNA™ Whole Genome Amplification Kits offer superior coverage and uniformity in single-cell WGS.

Ten single cells were isolated from a human B-lymphocyte cell culture (CEPH1463/NA12878/ GM12878 human genome reference standard). WGA was performed on individual cells, using the ResolveDNA Whole Genome Amplification Kit. WGA products were converted to indexed libraries and subjected to high-coverage whole genome sequencing (WGS) on the Illumina® platform. For the other WGA methods, low-pass WGS data (generated from individual BJ1 fibroblasts), were obtained from a previously published study (Chen, C. *et al. Science* 2017; 356: 189).

To achieve a fair comparison of the various WGA methods, raw data for all samples were aligned and pre-processed for variant calling using the same pipeline. All metrics shown in the table were generated from randomly subsampled BAM files (300 million reads per cell). Note that the metrics for all the methods other than the ResolveDNA kit are overestimates, due to the way in which data analysis was performed in the original study.

Bead Purification



The ResolveDNA™ Whole Genome Amplification and library preparation workflows include a number of bead-based purification steps. ResolveDNA Bead Purification Kits contain validated paramagnetic beads, as well as an optimized elution buffer to ensure highly efficient bead cleanups.

The ResolveDNA Dual Volume Strip Tube Magnet and Magnetic Plate are available to support all bead-based cleanups in the ResolveDNA™ Whole Genome Sequencing Workflow, in either plate or strip tube-based format.



The ResolveDNA Dual Volume Strip Tube Magnet is compatible with 8 strip 0.2 mL PCR tubes. One side handles volumes of 50 μ L – 0.2 mL and the other side 5 μ L – 50 μ L.




The ResolveDNA Magnetic Plate is compatible with full, semi and non-skirted PCR plates, and 0.2 mL PCR strips.

[Product Details >](#)

BioSkryb
GENOMICS

Library Preparation

OVERVIEW





To obtain superior whole genome sequencing data from single cells and ultra-low inputs, it is critical to preserve the quality of WGA reaction products generated with ResolveDNA™ WGA Kits during library preparation.

ResolveDNA Library Preparation Kits employ proprietary technology and reagent formulations to retain molecular diversity and introduce minimal bias during the construction of sequencing-ready libraries. The streamlined, optimized, ligation-based workflow does not require fragmentation of input DNA (WGA reaction products). ResolveDNA Multi-Use Library Adapters are supplied in a convenient plate-based format. These full-length adapters provide unique dual indices compatible with all Illumina(R) sequencers.

[Product Details >](#)

BioSkryb
GENOMICS

Sequencing



Appropriate experimental design is critical in single-cell genomics projects. Requirements may differ widely, depending on the application and study objectives.

The total amount of sequencing required for a study is determined by the size of the cell population to be interrogated, as well as the desired sequencing depth. BioSkryb offers high-quality sequencing services to ensure that the valuable information contained in every cell is preserved every step of the way; from sample preparation to analysis. Our whole genome sequencing service includes an initial (optional) round of low-pass sequencing to assess library quality to ensure the most robust libraries move forward to deep sequencing.

Contact our experienced team for assistance with experimental design, highly competitive sequencing services, or customized projects, including whole genome amplification, library preparation, sequencing and/or analysis. We are here to help you explore, discover and transform.

[Contact us >](#)

BioSkryb
GENOMICS

Analysis

OVERVIEW



BaseJumper

In addition to premium products and services, BioSkryb also offers innovative computational tools for analyzing human cells.

The cloud-based BaseJumper Bioinformatics Platform was developed to significantly shorten the time between data generation and insights.

The platform uses standardized best practices for sequencing analysis, and enables biologists (with no bioinformatics expertise) to explore the unique features of human single-cell data. Powerful visualization tools provide easy access to variant data, annotations, and quality metrics across multiple cells. Aligned data may also be downloaded for exploration using other applications.

**Contact us to learn more about
BaseJumper Bioinformatics Platform.**




BioSkryb
GENOMICS

Analysis

VISUALIZATION TOOLS




Circos plot and IGV viewer show variant density across the genome. Variants can be filtered by several criteria, including origin (somatic vs germline), synonymous vs non-synonymous, or the presence of COSMIC annotations.



Lineage browser allows the exploration of variants that likely occurred ancestrally. Branches may be highlighted and explored at higher resolution, via built-in annotation and link-outs to additional annotation sources.

Analysis

VISUALIZATION TOOLS



2D Cell representation shows similarity of cells based on their somatic variants.



Quality metrics enables easy visualization of QC metrics for a data set.

Ordering Information




STARTER PACKS & EQUIPMENT

Please contact us to request a quote or to speak to a member of our team. info@bioskryb.com

	PRODUCT CODE	PRODUCT DETAILS	DESCRIPTION	PACK SIZE
	100137	ResolveDNA™ Complete Starter Pack:	Everything needed to run PTA from start to finish. Each high-quality product in the ResolveDNA Complete Starter Pack has been carefully developed to provide optimal performance for the ResolveDNA Whole Genome Amplification Kit.	Bundle
		ResolveDNA™ Whole Genome Amplification Kit		96 reactions
		ResolveDNA™ Library Preparation Kit		96 reactions
		ResolveDNA™ Multi-Use Library Adapters		960 Reactions
		BioSkryb PCR Plate Spinner		1
		BioSkryb PCR Plate Thermal Mixer		1
		BioSkryb PCR Cooler		1
		BioSkryb Low Bind 96-well PCR Plates		Pack of 25
		PCR Plate Sealing Film		Pack of 100
		ResolveDNA™ PTA-Grade Cell Buffer Pack		12 x 550 µL
	100180	ResolveDNA™ Consumables Only Starter Pack:	The ResolveDNA Consumables Only Starter Pack includes consumables needed for optimal results with your ResolveDNA Whole Genome Amplification Kit.	Bundle
		ResolveDNA™ Whole Genome Amplification Kit		96 reactions
		ResolveDNA™ Bead Purification Kit		96 reactions
		ResolveDNA™ Library Preparation Kit		96 reactions
		ResolveDNA™ Multi-Use Library Adapters		960 reactions
		BioSkryb Low Bind 96-well PCR Plates		Pack of 25
		ResolveDNA™ PTA-Grade Cell Buffer Pack		12 x 550 µL
		PCR Plate Sealing Film		Pack of 100



Ordering Information

STARTER PACKS & EQUIPMENT

	PRODUCT CODE	PRODUCT DETAILS	DESCRIPTION	PACK SIZE
	100199	ResolveDNA™ FACS Kit: ResolveDNA™ PTA-Grade Cell Buffer BioSkryb Low Bind 96-well PCR Plates PCR Plate Sealing Film	The ResolveDNA FACS Kit contains 96-well plates and sealing film, as well as ResolveDNA PTA-Grade Cell Buffer—specifically optimized for PTA-based whole genome amplification from sorted cells.	Bundle 6 mL Pack of 25 Pack of 100
	100150	BioSkryb PCR Plate Thermal Mixer	The BioSkryb PCR Plate Thermal Mixer, controlled by a DC brushless motor and microcomputer, perfectly combines the constant temperature and shaking required during the PTA process, optimizing time and improving efficiency.	Each
	100153	BioSkryb PCR Plate Spinner	A high-quality plate spinner is required for quick spins throughout the PTA process. Designed for small samples in 96 or 384-well skirted, non-skirted and semi-skirted style plates.	Each





Ordering Information

WHOLE GENOME AMPLIFICATION KITS

	PRODUCT CODE	PRODUCT DETAILS	DESCRIPTION	PACK SIZE
	100068	ResolveDNA™ Whole Genome Amplification Kit	24 reaction PTA-based kit for accurate and reproducible whole genome amplification from single cells and ultra-low-input DNA inputs (4 pg to <10 ng).	24 reactions
	100136	ResolveDNA™ Whole Genome Amplification Kit	96 reaction PTA-based kit for accurate and reproducible whole genome amplification from single cells and ultra-low-input DNA inputs (4 pg to <10 ng).	96 reactions

Ordering Information



BEAD PURIFICATION KITS AND MAGNETS

	PRODUCT CODE	PRODUCT DETAILS	DESCRIPTION	PACK SIZE
	100121	ResolveDNA™ Bead Purification Kit	For reliable cleanup of WGA reaction products and libraries prepared for Illumina® sequencing.	24 reactions
	100182	ResolveDNA™ Bead Purification Kit	For reliable cleanup of WGA reaction products and libraries prepared for Illumina® sequencing.	96 reactions
	100135	ResolveDNA™ Magnetic Plate*	Compatible with full-skirted, semi-skirted and non-skirted PCR plates and 0.2 mL PCR strip tubes.	Each
	100226	ResolveDNA™ Dual Volume Strip Tube Magnet	Compatible with 8-strip 0.2 mL PCR tubes. Use one side of the magnet for volumes from 5 µL – 50 µL, and the other side for 50 µL – 0.2 mL volumes.	Each

Ordering Information

LIBRARY PREPARATION REAGENTS

Please contact us to request a quote or to speak to a member of our team. info@bioskryb.com

	PRODUCT CODE	PRODUCT DETAILS	DESCRIPTION	PACK SIZE
	100080	ResolveDNA™ Library Preparation Kit	Optimized ligation-based (without fragmentation) conversion of PTA products to Illumina®-ready libraries for sequencing.	96 reactions
	100181	ResolveDNA™ Multi-Use Library Adapters	Unique dual-indexed adapters compatible with all Illumina® sequencers. Compatible with the ResolveDNA™ Whole Genome Amplification Kit at a 10x dilution, allowing for multiple uses per adapter plate.	960 Reactions

bioskryb.com

All data on file.
For Research Use Only. Not for use in diagnostic procedures.
BIOSKRYB, RESOLVEDNA and BASEJUMPER are trademarks of BioSkryb, Inc.
All other product names and trademarks are the property of their respective owners.

© 2020. BioSkryb, Inc. All rights reserved.

