

Application: Cell and Gene Therapy (CGT)



Empower Your CGT Development Programs with Sensitive and Precise Single-Cell Multiomics

Challenges in CGT Development



Viral Gene Therapies

- Insertional mutagenesis can give rise to neoplastic changes
- Uncertainty in transduction efficiency, gene expression, and tissue targeting



CRISPR/Cas9 Gene Editing

- Off-target Cas9 activity results in mutations that vary between cells
- Difficult to confirm zygosity




Induced Pluripotent Stem Cell (iPSC) Therapies

- Oncogenic potential from genomic and transcriptomic variation in iPSCs
- Heterogeneity in cell identity


Solutions from BioSkryb Genomics

Primary template-directed amplification (PTA) technology, available in **ResolveDNA[®]** and **ResolveOME[™]** kits, enables multiomic analysis at the single-cell level, empowering CGT researchers to sensitively and precisely detect both desirable and deleterious changes in individual candidate cells.



Detect off-target events in single cells comprehensively.

Analysis of PTA-enabled whole genome sequencing data¹ reveals hundreds of indels detected following CRISPR/Cas9 gene editing with EMX1 or VEGFA guide RNAs in single cells (SC). Total number of indels detected versus predicted in each group shown.



Identify underlying mechanisms of gene expression changes in individual cells.

Analysis of parental and drug-resistant cell lines following ResolveOME whole genome and transcriptome amplification reveals single nucleotide variants (SNVs) in non-coding upstream regulatory regions that are associated with changes in gene expression.²

References:

1. Gonzalez-Pena et al., "Accurate genomic variant detection in single cells with primary template-directed amplification". PNAS. (2021).
2. Marks et al., "Unifying comprehensive genomics and transcriptomics in individual cells to illuminate oncogenic and drug resistance mechanisms". bioRxiv. (2023)

Assay Performance

Table 1: ResolveOME WGS DNA Performance

Characteristic	Observed Values
Accuracy	99.99%
Sensitivity	96.65%
Specificity	99.99%
Allelic Balance	91.20%
Genomic Coverage	97.59%


Table 2: ResolveOME WTS RNA Performance

Characteristic	Observed Values
Protein Coding Genes	3451 ± 732
Concordance	0.97
Variance (CV)	32.9%

Assay performance characteristics of DNA (Table 1) and RNA (Table 2) isolated using ResolveOME Whole Genome and Transcriptome Single-Cell Core Kit. Analysis of FACS-sorted NA12878 single cells prepared with ResolveOME versus gold-standard reference. WGS: whole genome sequencing. WTS: whole transcriptome sequencing.

ResolveServicesSM

Custom-built service projects, from singulating cells to figures. Services can include:



FACS: Fluorescence-Activated Cell Sorting, WGA: Whole Genome Amplification, WTA: Whole Transcriptome Amplification, **Optional **targeted protein** detection available

Products

Product	Description	Codes	Well Format
ResolveOME™ Whole Genome and Transcriptome Single-Cell Core Kit	PTA-based kit for whole genome and transcriptome amplification plus NGS library preparation from single cells.	100956	96
		100957	384
ResolveDNA® Whole Genome Single-Cell Core Kit	PTA-based kit for whole genome amplification plus NGS library preparation from single cells.	100954	96
		100955	384
ResolveDNA® Whole Genome Amplification Kit	PTA-based kit for whole genome amplification from as little as 4 pg input DNA.	100545	Up to 384*
BaseJumper® Bioinformatics Platform	A platform for multiomic data analysis.	100605	-
ResolveXOME™ Exome Capture Module	Exome capture module for use with ResolveDNA or ResolveOME Single-Cell Core Kits.	Early-Access	-

*Dependent on workflow option used

For a complete list of services, products, and pricing, email a member of our team, info@bioskryb.com



All data on file.
 BIOSKRYB, RESOLVEDNA, RESOLVEOME, RESOLVEOME, RESOLVEOME, RESOLVESERVICES, and BASEJUMPER are trademarks of BioSkryb, Inc.
 All other product names and trademarks are the property of their respective owners.
 © 2024 BioSkryb, Inc. All Rights Reserved

