

# Application: Hematologic Malignancies

## Transform Your Hematology Research with Sensitive and Precise Single-Cell Multiomics

### Challenges in Hematologic Oncology Research



#### Rare malignant clones and persister cells facilitate resistance and recurrence

Bulk sequencing may miss significant mutations in rare cells and fail to identify which mutations co-occur in the same cell.



#### Hematologic malignancies are driven by various sizes of genomic alterations

Panel-based DNA approaches struggle to capture single nucleotide variants (SNV), copy number variation (CNV), and structural variants (SV) in one assay.




#### Consequences of SNV, CNV, and SV changes are not always obvious

DNA-only approaches can fail to identify the gene expression changes resulting from mutations in both malignant and benign hematologic disorders.

### Solutions from BioSkrbyb Genomics

**ResolveDNA®** enables single-cell whole genome sequencing, empowering researchers to sensitively and precisely detect whole genome CNV and SNVs at a single cell resolution. **ResolveOME™** allows investigators to move beyond inference from gene expression and define the genetic mechanisms controlling gene expression at a single-cell level.



**Detect SNV changes in individual therapy-resistant cells.** Analysis of parental and quizartinib-resistant cell lines following ResolveOME whole genome and transcriptome amplification. **A.** Previously published *FLT3* internal tandem duplications and N841K mutations are detected. **B.** SNVs in individual cells vary between sensitive and resistant lines.<sup>1</sup>

**Identify underlying mechanisms of gene expression changes in individual cells.** Analysis of parental and drug-resistant cell lines following ResolveOME whole genome and transcriptome amplification reveals SNVs in non-coding upstream regulatory regions that are associated with changes in gene expression.<sup>1</sup>

### Reference

1. Marks et al., "Unifying comprehensive genomics and transcriptomics in individual cells to illuminate oncogenic and drug resistance mechanisms". *bioRxiv*. (2023)

For Research Use Only. Not for Use in Diagnostic Procedures.

## ResolveDNA® and ResolveOME™ Assay Performance

**Table 1: ResolveOME WGS DNA Performance\***

Characteristic	Observed Values
Accuracy	99.5%
Sensitivity	97.1%
Specificity	99.2%
Allelic Balance	98.4%
Genomic Coverage	97.1%

**Table 1: Assay performance characteristics of DNA isolated using ResolveOME.** Analysis of FACS-sorted NA12878 single cells prepared with ResolveOME versus gold-standard reference. WGS: whole genome sequencing.

\*DNA amplified using ResolveDNA and ResolveOME have comparable DNA performance characteristics. All data on file.

**Table 2: ResolveOME WTS RNA Performance**

Characteristic	Observed Values
Genes Detected	4,546
Reportable Range	6,057
Average Concordance	0.91
Reproducibility (CV)	43.3%

**Table 2: Assay performance characteristics of RNA isolated using ResolveOME.** Analysis of FACS-sorted NA12878 single cells prepared with ResolveOME versus gold-standard reference. WTS: whole transcriptome sequencing.

### Custom Services

We offer custom service packages from our end-to-end single-cell multiomic pipeline, from singulating cells to ready-to-publish figures. All services include quality control verification. Services include:

- Cell sorting from fresh or frozen cells and tissues
- Whole genome amplification or whole genome and transcriptome amplification
- Library preparation for downstream applications, such as whole genome or targeted sequencing
- Sequencing of 550M quality reads capturing >97% of the human genome from each cell
- Analysis using our bioinformatics platform, **BaseJumper™**

### Products

Code	Product	Description
100500	ResolveOME™ Whole Genome and Transcriptome Amplification System	PTA-based kit for accurate and reproducible whole genome and transcriptome amplification from single cells.
100545	ResolveDNA® Whole Genome Amplification Kit	PTA-based kit for accurate and reproducible whole genome amplification from single cells and low-input DNA inputs.
100605	BaseJumper™ Bioinformatics Platform	A complete bioinformatics solution for multiomic data analysis and visualization.

*For a complete list of services, products, and pricing, email a member of our team, [info@bioskryb.com](mailto:info@bioskryb.com)*

**BioSkryb**  
GENOMICS

All data on file.

BIOSKRYB, RESOLVEDNA, RESOLVEOME, and BASEJUMPER are trademarks of BioSkryb, Inc.

All other product names and trademarks are the property of their respective owners.

© 2023 BioSkryb, Inc. All Rights Reserved

TAS\_059 | 11/2023



[bioskryb.com/resources](https://bioskryb.com/resources)

For Research Use Only. Not for Use in Diagnostic Procedures.