


I. Salas-González<sup>1</sup>, T. Morozova<sup>1</sup>, D. Arvapalli<sup>1</sup>, S. Velivela<sup>1</sup>, G. Harton<sup>1</sup>, J. Zawistowski<sup>1</sup>, J. Blackinton<sup>1</sup>, V. Weigman<sup>1</sup>, J. A. A. West<sup>1</sup>  
<sup>1</sup>BioSkyrb Genomics, Inc., Durham, NC

## Abstract

Single-cell RNAseq methodologies have revolutionized the ability to define cell type identity and transcriptional phenotypic diversity within a tumor sample, yet in most studies DNA-level variation contributing to driving this diversity of gene expression remains uncharacterized. We employed ResolveOME chemistry to bridge this gap within an individual cell, by taking inventory of genome-wide single nucleotide variation (SNV) and cross referencing it with full-transcript RNAseq data in an acute myeloid leukemia (AML) cell line model of drug resistance. The catalog of genome-wide SNV in conjunction with transcriptomic data in this model powered an association screening to identify regulatory SNVs with biased prevalence in single cells resistant to the FLT3 inhibitor quizartinib vs treatment-naïve single cells, which then were correlated with differential expression of genes proximal to the variant in the same cell. Initially, we focused on nucleotide variation within core promoter regions and within gene bodies. Fitting a zero-inflated linear model to a matrix of expression and genotype across all single cells revealed an intronic heterozygous variant in MYC in parental single cells, while, in contrast, the alternate allele was absent in quizartinib-resistant cells. Upregulation of MYC transcript was correlated with the absence of this variant in resistant single cells, suggesting that the single nucleotide change may have intronic enhancer activity as opposed to intronic regulation of splicing. Additionally, the screen uncovered a heterozygous single nucleotide change within 5kb 5' of the transcriptional start site of the mRNA binding factor PABPC4 as a candidate promoter variant not present in parental single cells yet harbored by 50% of the quizartinib resistant single cells displaying PABPC4 expression relative to parental cells. Differential expression analysis uncovered the enhancer factor CEBPA as upregulated in the resistant single cells, and the genomic component of ResolveOME identified two SNVs 20 kb upstream of the locus representing putative CEBPA enhancer variants influencing the differential gene expression. The co-identification of CEBPA and PABPC4 through this multi-omic approach suggest that quizartinib resistance in this model is mediated in part by global gene regulation through enhancer modulation and through mRNA stability/translational regulation, respectively. We are furthering enhancer variant detection in intergenic space by overlaying ChIP-seq, chromatin accessibility, and transcription factor binding site data to prioritize candidates contributing to resistance in this AML drug resistance model.


## Drug Resistance: linking variants' impact on expression

To demonstrate the utility of concurrent genomic and transcriptomic information in single cells in the context of drug resistance, we created a model by exploiting the presence of an internal tandem duplication (ITD) mutation in MOLM-13 cells. Since the ITD mutation, found in ~20% of AML patients, hyperactivates FLT3 signaling and results in poor prognosis and relapse, we treated non-resistant (Parental) and drug-sensitive cells with a continual dose of 2 nM quizartinib (Resistant).




## Drug Resistance: expression impacted by gene copy

Given the concurrent characterization of DNA and RNA from the same cell, we had the power to associate genome-wide the influence of copy number (CNV) events with respect to gene expression. Interestingly, we identified clusters of genes that showed clear positive and negative association with ploidy across the dataset. Further, we also identified genes with more cryptic patterns of influence suggesting alternative mechanisms influencing gene expression such as epigenetics.




## Tumor Heterogeneity: characterizing cell states

We demonstrate analogous multi-omic utility in elucidating single-cell oncogenic mechanisms in primary human cancer. To do so, we characterized genomic and transcriptomic contributions to the transition of premalignant carcinoma in situ (DCIS) to invasive ductal carcinoma (IDC).



## Drug Resistance: Differential isoform presentation

The ability of ResolveOME™ to characterize the full-length transcriptome of a single cell opens the possibility to identify differential transcript usage (DTU) events that could be driving physiological phenotypes. To do so, we performed DTU characterization in the MOLM-13 model described previously. Briefly, a DTU event across two treatment groups (e.g. Parental and Resistant MOLM-13 cells) tests for proportional differences in the expressed transcript composition of a gene, thus comparing how much each transcript contributes to the gene's total expression between conditions.



## ResolveOME workflow

ResolveOME permits to simultaneously characterize the genomic and transcriptomic data from the same individual cell vastly increases the complexity of putative mechanisms of drug resistance and oncogenesis

