

Application: Cell and Gene Therapy (CGT)

Empower Your CGT Development Programs with Sensitive and Precise Single-Cell Multiomics

Challenges in CGT Development



Viral Gene Therapies

- Insertional mutagenesis can give rise to neoplastic changes
- Uncertainty in transduction efficiency, gene expression, and tissue targeting



CRISPR/Cas9 Gene Editing

- Off-target Cas9 activity results in mutations that vary between cells
- Difficult to confirm zygosity




Induced Pluripotent Stem Cell (iPSC) Therapies

- Oncogenic potential from genomic and transcriptomic variation in iPSCs
- Heterogeneity in cell identity

Solutions from BioSkryb Genomics

Primary template-directed amplification (PTA) technology, available in [ResolveDNA®](#) and [ResolveOME™](#) kits, enables multiomic analysis at the single-cell level, empowering CGT researchers to sensitively and precisely detect both desirable and deleterious changes in individual candidate cells.




Detect off-target events in single cells comprehensively.

Analysis of PTA-enabled whole genome sequencing data¹ reveals thousands of indels detected following CRISPR/Cas9 gene editing with EMX1 or VEGFA guide RNAs in single cells (SC). Proportion of or total number of indels detected versus predicted in each group shown.

References

1. Gonzalez-Pena et al., "Accurate genomic variant detection in single cells with primary template-directed amplification". *PNAS*. (2021).
2. Marks et al., "Unifying comprehensive genomics and transcriptomics in individual cells to illuminate oncogenic and drug resistance mechanisms". *bioRxiv*. (2023)



Identify underlying mechanisms of gene expression changes in individual cells.

Analysis of parental and drug-resistant cell lines following ResolveOME whole genome and transcriptome amplification reveals single nucleotide variants (SNVs) in non-coding upstream regulatory regions that are associated with changes in gene expression.²

ResolveDNA® and ResolveOME™ Assay Performance

Table 1: ResolveOME WGS DNA Performance*

Characteristic	Observed Values
Accuracy	99.5%
Sensitivity	97.1%
Specificity	99.2%
Allelic Balance	98.4%
Genomic Coverage	97.1%

Table 2: ResolveOME WTS RNA Performance

Characteristic	Observed Values
Genes Detected	4,546
Reportable Range	6,057
Average Concordance	0.91
Reproducibility (CV)	43.3%

Table 1: Assay performance characteristics of DNA isolated using ResolveOME. Analysis of FACS-sorted NA12878 single cells prepared with ResolveOME versus gold-standard reference. WGS: whole genome sequencing.

Table 2: Assay performance characteristics of RNA isolated using ResolveOME. Analysis of FACS-sorted NA12878 single cells prepared with ResolveOME versus gold-standard reference. WTS: whole transcriptome sequencing.

*DNA amplified using ResolveDNA and ResolveOME have comparable DNA performance characteristics. All data on file.

Custom Services

We offer custom service packages from our end-to-end single-cell multiomic pipeline, from singulating cells to ready-to-publish figures. All services include quality control verification. Services include:

- Cell sorting from fresh or frozen cells and tissues
- Whole genome amplification or unified whole genome and transcriptome amplification
- Library preparation for downstream applications, such as whole genome or targeted sequencing
- Sequencing of 550M quality reads capturing >97% of the human genome from each cell
- Analysis using our bioinformatics platform, **BaseJumper™**

Products

Code	Product	Description
100500	ResolveOME™ Whole Genome and Transcriptome Amplification System	PTA-based kit for accurate and reproducible whole genome and transcriptome amplification from single cells.
100545	ResolveDNA® Whole Genome Amplification Kit	PTA-based kit for accurate and reproducible whole genome amplification from single cells and low-input DNA inputs.
100605	BaseJumper™ Bioinformatics Platform	A complete bioinformatics solution for multiomic data analysis and visualization.

For a complete list of services, products, and pricing, email a member of our team at info@bioskryb.com

BioSkryb
GENOMICS

All data on file.

BIOSKRYB, RESOLVEDNA and BASEJUMPER are trademarks of BioSkryb, Inc.

All other product names and trademarks are the property of their respective owners.

© 2023 BioSkryb, Inc. All Rights Reserved

TAS_055 | 09/2023



bioskryb.com

For Research Use Only. Not for Use in Diagnostic Procedures.