

Unlock Your Discovery Potential with Single-Cell Whole Genome Amplification

ResolveDNA Whole Genome Amplification System



Enables whole genome, exome, or panel sequencing from a single cell

Whole genome amplification from single cells or low-input (4 pg to <10 ng) DNA samples



Provides industry-leading genomic coverage

Leverages a novel patented technology, primary template-directed amplification (PTA)¹, to dramatically increase genomic capture and coverage to 97%



Results in unprecedented single nucleotide variant (SNV) calling resolution

Reduces biases, experimental artifacts, and poor reproducibility associated with other whole genome amplification methods²



Fits into established laboratory protocols

Compatible with various methods of single-cell singulation, sequencing platforms, and downstream applications, including whole exome and panel-based sequencing³



Scales to experiment size

Low-cost, automation-compatible, scalable approach with up to 384 reactions per kit



Bioinformatics analysis available with BaseJumper[®]

Bioinformatics analysis and data visualization through BaseJumper bioinformatics platform and computational biology custom services



References:

1. Gonzalez-Pena V, et al. Proc. Natl. Acad. Sci. U.S.A. 2021; 118 (24): e2024176118; doi: 10.1073/pnas.2024176118
2. Luquette L, et al. Nat Gen. 2022; 54: 1564–1571. doi: 10.1038/s41588-022-01180-2
3. Data on file

Assay Performance

ResolveDNA WGS Performance	
Characteristic	Observed Values
Accuracy	99.99%
Sensitivity	96.90%
Specificity	99.99%
Allelic Balance	93.80%
Genomic Coverage	97.50%

Figure 1: Assay Performance Characteristics.

Analysis of NA12878 single cells prepared with ResolveDNA Whole Genome Single-Cell Core Kit versus gold standard reference sequence. WGS: whole genome sequencing.

A Revolution in Resolution From Each Cell


DNA

- Resolve SNV
- Resolve SV
- Resolve CNV
- Resolve Ploidy
- Resolve Exomes
- Resolve Genomes
- Resolve Panels
- Resolve Edits

Resolve More.

ResolveServicesSM

Custom-built service projects, from singulating cells to figures. Services can include:



FACS: Fluorescence-Activated Cell Sorting, WGA: Whole Genome Amplification, WTA: Whole Transcriptome Amplification, **Optional targeted protein detection available

Products

Product	Description	Codes	Well Format
ResolveOME™ Whole Genome and Transcriptome Single-Cell Core Kit	PTA-based kit for whole genome and transcriptome amplification plus NGS library preparation from single cells.	100956	96
		100957	384
ResolveDNA® Whole Genome Single-Cell Core Kit	PTA-based kit for whole genome amplification plus NGS library preparation from single cells.	100954	96
		100955	384
ResolveDNA® Whole Genome Amplification Kit	PTA-based kit for whole genome amplification from as little as 4 pg input DNA.	100545	Up to 384*
BaseJumper® Bioinformatics Platform	A platform for multiomic data analysis.	100605	-
ResolveXOME™ Exome Capture Module	Exome capture module for use with ResolveDNA or ResolveOME Single-Cell Core Kits.	Early-Access	-

*Dependent on workflow option used

For a complete list of services, products, and pricing, email a member of our team, info@bioskryb.com

BioSkryb
GENOMICS

bioskryb.com

All data on file.
BIOSKRYB, RESOLVEDNA, RESOLVEOME, RESOLVEXOME, RESOLVESERVICES, and BASEJUMPER are trademarks of BioSkryb, Inc.

All other product names and trademarks are the property of their respective owners.

© 2025 BioSkryb, Inc. All Rights Reserved

TAS_044.6 | 02/2025

For Research Use Only. Not for Use in Diagnostic Procedures.

